
Saline-Air Device
Dispositif à air et sérum physiologique
Dispositivo a soluzione salina-aria
Kochsalzlösung-Luft-Gerät
Dispositivo de solución salina-aire
Saltvand-luft anordning
Zoutoplossing-luchtapparaat
Anordning för koksaltlösning-luft
Salttidig luft-apparat
Dispositivo de ar-solução salina


Instructions for Use
Mode d'emploi
Istruzioni per l'uso
Gebrauchsanweisung
Instrucciones de uso
Brugsanvisning
Gebruiksaanwijzing
Bruksanvisning
Bruksanvisning
Instruções de uso

English for EU only	3
Français	6
Italiano.....	9
Deutsch	12
Español.....	15
Dansk.....	18
Nederlands	21
Svenska.....	24
Norsk	27
Português	30

Instructions for Use

DEVICE DESCRIPTION

The FemVue Saline-Air Device (FemVue) is a dual-barrel contrast media syringe that can be connected to an intrauterine catheter to instill saline-air contrast media during sono-hysterosalpingogram (Sono HSG) procedures. Sono HSG consists of an ultrasound evaluation of the fallopian tubes with or without assessment of the uterine cavity.

INDICATIONS FOR USE

The FemVue Saline-Air Device instills a consistent alternating pattern of saline and air as a continuous stream of contrast media into the uterus and fallopian tubes to be used in conjunction with an intrauterine catheter for performance of sono-hysterosalpingogram (Sono HSG).

CONTRAINDICATIONS

The FemVue Saline-Air Device should not be used in any woman who has a contraindication to hysterosalpingography, including women who are pregnant or who have been pregnant in the previous 6 weeks (including miscarriage). These women may be at increased risk for air embolism.

HOW SUPPLIED

Sterile for single use only.

STORAGE

Store in a cool, dry place.

WARNINGS/PRECAUTIONS

- Do not use oil-based contrast media. Saline-air is the recommended contrast media for this device.
- In order to minimize the risk of air embolism, do not exceed delivery of six (6) filled syringe volumes to the patient. Air embolism has not been reported with saline-air contrast volumes below 70 mL.
- Sono HSG with FemVue should be performed after completion of the menstrual cycle and before the onset of ovulation.
- Intended for single patient use only. Reuse creates a potential risk of patient or user infections.
- Do not use if pouch is damaged.
- When performing Sono HSG in a patient who has had a tubal occlusion procedure for permanent contraception, refer to the tubal occlusive device's Instructions for Use for confirmation testing.

CLINICAL USE

The table below presents reported sensitivity, specificity, and concordance of saline-air sono-hysterosalpingography to assess tubal patency compared to laparoscopic chromopertubation. The FemVue Saline-Air Device has not been evaluated in a clinical study.

Reference	Number of Tubes	Sensitivity (%)	Specificity (%)	Concordance (%)
Allahbadia Fertility and Sterility. 1992;58(5):901-907.	134	Not reported		94
Heikkinen et al. Fertility and Sterility. 1995;64(2):293-298.	61	82	86	85
Volpi et al. Ultrasound Obstetrics Gynecology. 1996;7:43-48.	56	85	92	89
Spalding et al. Human Reproduction. 1997;12(11):2461-2464.	29	63	100	90
Inki et al. Acta Obstetrica et Gynecologica Scandinavica. 1998;77:978-982.	53	90	83	89
Jeanfy et al. Journal Ultrasound Medicine. 2000;19:519-527.	29	86	77	79
Exacoustos et al. The Journal of the American Association of Gynecologic Laparoscopists. 2003;10(3):367-372.	30	75	91	87

FEMVUE PREPARATION

1. Completely submerge tip into a saline-filled bowl.
2. Fully pull back the plunger handle and hold until the saline chamber completely fills. Do not squeeze the air and saline chambers.

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

1. Sonographic Landmarking.

Conduct your standard ultrasound evaluation per your practice guidelines.

Attempt to locate the following in the **transverse view**:

- Endometrial stripe and uterotubal junctions (**Fig A**)
- Position of each ovary relative to the uterus (**Fig B**)

Take note of the probe position.


When performing Sono HSG in a patient who has had a tubal occlusion procedure for permanent contraception, first follow any applicable Instructions for Use of the tubal occlusive device.

2. Insert Intrauterine Catheter per the Catheter Instructions for Use.

If desired, first perform Saline Infusion Sonography with a saline-filled syringe per your practice protocol, because bubbles in the uterus from the FemVue may cause artifact.

3. Inflate Balloon, if applicable.

Inflation of the balloon with subsequent placement over the internal cervical os is recommended to prevent retrograde flow. (**Fig C**)


4. Prime FemVue and Connect to Catheter. (**Fig D**)

To prime FemVue, submerge the tip in the saline-filled bowl and depress the plunger handle until bubbles are visible.

Ensure FemVue is primed just before attachment to avoid delay in contrast visualization.

Do not overtighten FemVue's luer to the catheter. This prevents catheter kinking and luer detachment while ensuring easy device removal for refilling, if necessary.

5. Deliver Contrast Slowly.

With ultrasound probe in place, slowly and steadily depress the plunger handle to deliver the contrast. Visualize saline-air contrast (bubbles) entering the cavity (distention is not necessary).

Confirm in the **sagittal view** there is no retrograde flow through the cervix. If needed, adjust the balloon's placement or use a balloon to block the flow. (**Fig E**)

Start with one fill of the FemVue Saline-Air Device. Use the minimum number of fills necessary to perform tubal assessment. Do not exceed six (6) filled syringe volumes.

Troubleshooting: No contrast exiting catheter and plunger resistance

- Make certain that catheter clamp is open, if applicable.
- Confirm the catheter is patent.
- Ensure FemVue has not been over tightened to catheter luer.
- Replace catheter if kinking is suspected.

6. Make Tubal Assessment.

In the **transverse view** orient the probe to observe the uterotubal junction to assess contrast flow in one tube. Although flow may be seen bilaterally, evaluate each tube individually. (**Fig F**)

- Locate flow in Zone 1 and hold view to observe.
- Slowly and methodically scan to possibly observe tubal flow in the remaining zones.

Evaluate contralateral tube.

Troubleshooting: No contrast flow visible in Zone 1

- Ensure there is no retrograde flow through the cervix.
- Consider repositioning ultrasound probe.
- Hold probe, maintain plunger handle position, and wait for resolution of possible tubal spasm.
- After extensive observation, trace laterally to Zone 3 and observe.

Troubleshooting: Inadequate visualization of tube

- Consider pulling back with probe if tube is viewed anteriorly.
- Consider adjusting ultrasound equipment settings.

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Guidelines from published literature using saline and air are as follows:

- **Zone 1** - Flow in the interstitial part of the tube (minimum criterion).^{1,2}
- **Zone 2** - Flow throughout course of tube (may not be seen).¹
- **Zone 3** - Flow exiting tube (fimbrial turbulence, bubbles seen around ovary or in cul-de-sac).^{1,2,3}



SONO HSG: DEMONSTRATING TUBAL OCCLUSION WITH PERMANENT CONTRACEPTIVE DEVICE IN PLACE

Guidelines from published literature using sonographic contrast are as follows:

- **Zone 1** - Flow to occlusive device, and
- **Zone 2** - Absence of flow through and beyond occlusive device, and
- **Zone 3** - Absence of flow exiting tube (fimbrial turbulence, bubbles around ovary and/or in cul-de-sac).^{4,5,6}

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© 2012 Femasys Inc.
Femasys Inc.